

Datasheet

Nitrotyrosine polyclonal antibody (HRP)

Catalog Number: PAB14356

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against Nitrotyrosine.

Immunogen: Nitrotyrosine conjugated with KLH.

Host: Goat

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody can be used for detection of nitrotyrosine in plasma and lipoproteins.

Form: Liquid

Conjugation: HRP

Recommend Usage: ELISA (1:100-1:4000)

Western Blot (1:100-1:4000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM PBS, 100 mM NaCl, pH 7.4
(10 mg/mL BSA, 0.01% thimerosal)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.