

Datasheet

RYBP purified MaxPab mouse polyclonal antibody (B02P)

Catalog Number: H00023429-B02P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human RYBP protein.

Immunogen: RYBP (NP_036366.3, 1 a.a. ~ 228 a.a) full-length human protein.

Sequence:

MTMGDKKSPTRPKRQAKPAADEGFWDCSVCTFRNSA
EAFKCSICDVRKGTSTRKPRINSQLVAQQVAQQYATP
PPPKEKKEKVEKQDKEKPEKDKEISPSVTKNNTNKK
KPKSDILKDPPSEANSIQSANATTKTSETNHTSRPRLK
NVDRSTAQQLAVTVGNVTVIITDFKEKTRSSSTSSSTV
TSSAGSEQNQSSSGSESTDKGSSRSSTPKGDMSAV
NDESF

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 23429

Gene Symbol: RYBP

Gene Alias: AAP1, DEDAF, YEAF1