

Datasheet

Collagen Type V polyclonal antibody (Biotin)

Catalog Number: PAB14400

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against native COL5.

Immunogen: Native human COL5.

Host: Goat

Reactivity: Human

Applications: Dot, ELISA, ICC

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Reacts with conformational determinants on human type V collagen as demonstrated by ELISA. May react with type V collagen from other species. Exhibits <10% cross reactivity with collagen type I, II, III, IV and VI. This antibody has not been tested for reactivity with other ECM proteins (e.g., laminin, fibronectin).

Form: Liquid

Conjugation: Biotin

Recommend Usage: Immunocytochemistry (1:20-1:40)

Immunohistochemistry (1:20-1:40)

ELISA (1:1000-1:4000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C.