

## Datasheet

### Human immunodeficiency virus P24 polyclonal antibody

**Catalog Number:** PAB14679

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against recombinant human Immunodeficiency Virus (HIV) P24.

**Immunogen:** Recombinant protein corresponding to human Immunodeficiency Virus (HIV) P24.

**Host:** Rabbit

**Reactivity:** Human, Viruses

**Applications:** ELISA

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody recognizes both native and recombinant HIV-P24.

**Form:** Lyophilized

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 20 mM PBS, pH 7.4

**Storage Instruction:** Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.