

Datasheet

Ulex Europaeus Lectin 1 polyclonal antibody

Catalog Number: PAB15537

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against native Ulex Europaeus Lectin 1.

Immunogen: Native purified from Ulex europaeus agglutinin isolated from seeds of Ulex europaeus (gorse seed).

Host: Rabbit

Reactivity: Human, Plants

Applications: IF, IHC, IHC-Fr, IHC-P
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody may be used for studies of the vascular structures. It reacts specifically with UEA-I bound to human endothelial cells of normal and neoplastic blood and lymphatic vessels. Its uniform reaction with vessels enables differentiation between vessels and artificial slits in the tissue. Other species have not been tested.

Form: Liquid

Recommend Usage: Immunofluorescence (1:500)
Immunohistochemistry (Frozen sections) (1:500)
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:500)
The optimal working dilution should be determined by the end user.

Storage Buffer: In antiserum (15 mM sodium azide)

Storage Instruction: Store at 4°C for short term. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.