

Datasheet

RAD51 polyclonal antibody

Catalog Number: PAB15573

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against full length recombinant RAD51.

Immunogen: Recombinant protein corresponding to full length *Saccharomyces cerevisiae* RAD51.

Host: Rabbit

Reactivity: Yeast

Applications: IF, IP, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (1:2000-1:10000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In antiserum (0.05% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 856831

Gene Symbol: RAD51

Gene Alias: MUT5