

Datasheet

Dnmt1 polyclonal antibody

Catalog Number: PAB15591

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against partial recombinant Dnmt1.

Immunogen: Recombinant protein corresponding to amino acids 1037-1386 mouse Dnmt1.

Host: Rabbit

Reactivity: Frog, Human, Mouse

Applications: IF, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (0.2-1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In antiserum (0.05% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 13433

Gene Symbol: Dnmt1

Gene Alias: Cxxc9, Dnmt, Dnmt1o, MTase, Met-1, Met1, MommeD2