

Datasheet

Hepatitis C virus NS-3 monoclonal antibody, clone 328

Catalog Number: MAB6333

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against recombinant hepatitis C virus NS-3.

Clone Name: 328

Immunogen: Recombinant protein corresponding to Human hepatitis C virus (HCV) NS-3.

Host: Mouse

Reactivity: Human, Viruses

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with recombinant NS-3 (amino acids 1252-1477) and recombinant Chimeric HCV Polyprotein (60 KDa).

Form: Lyophilized

Isotype: IgG1, kappa

Recommend Usage: Western Blot (0.1-1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM PBS, pH 7.2

Storage Instruction: Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.