

## Datasheet

### Hepatitis C virus NS-4 monoclonal antibody, clone 101

**Catalog Number:** MAB6342

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against recombinant hepatitis C virus NS-4.

**Clone Name:** 101

**Immunogen:** Recombinant protein corresponding to hepatitis C virus NS-4.

**Host:** Mouse

**Reactivity:** Human, Viruses

**Applications:** ELISA, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody reacts with synthetic NS-4a Protein (amino acids 1689-1735) and Recombinant Chimeric HCV Polyprotein (60 KDa).

**Form:** Lyophilized

**Isotype:** IgG1, kappa

**Recommend Usage:** Western Blot (0.5 ug/mL)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 10 mM PBS, pH 7.0

**Storage Instruction:** Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.