

Datasheet

Hepatitis E virus ORF3 monoclonal antibody, clone S-13

Catalog Number: MAB6345

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant hepatitis E virus ORF3.

Clone Name: S-13

Immunogen: Recombinant protein corresponding to least 30 amino acids in length of Chinese hepatitis E virus ORF3 strain.

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Isotype: IgM

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM PBS, pH 7.0

Storage Instruction: Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.