

## Datasheet

### Human immunodeficiency virus type 1 P24 (core protein) monoclonal antibody, clone 473

**Catalog Number:** MAB6346

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against recombinant human immunodeficiency virus type 1 P24 (core protein).

**Clone Name:** 473

**Immunogen:** Recombinant protein corresponding to human immunodeficiency virus type 1 P24 core protein.

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody reacts with HIV-1 P24 (core protein).

**Form:** Lyophilized

**Isotype:** IgG1, kappa

**Recommend Usage:** Western Blot (0.1-1.0 ug/mL)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 10 mM PBS, pH 7.2

**Storage Instruction:** Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.