

Datasheet

BXDC1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00084154-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human BXDC1 protein.

Immunogen: BXDC1 (NP_115570.1, 1 a.a. ~ 306 a.a) full-length human protein.

Sequence:

MDTLDRVVKPKTKRAKRFLEKREPKLNENIKNAMLIK
GNANATVTKVLKDVYALKKPYGVLYKKKNITRPFEDQT
SLEFFSKSDCSLFMFGSHNKKRPNNLVIGRMYDYHV
LDMIELGIENFVSLKDIKNSKCPEGTKPMLIFAGDDFDV
TEDYRRLKSLIDFFRGPTVSNIRLAGLEYVLHFTALNG
KIYFRSYKLLKKSGCRTPRIELEEMGPSLDLVLRRTHL
ASDDLKLSMKMPKALKPKKKKNISHDTFGTTYGRIHM
QKQDLSKLQTRKMKGLKKRPAERITEDHEKSKRIKKN

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84154

Gene Symbol: BXDC1

Gene Alias: FLJ21087, RPF2, bA397G5.4