

Datasheet

GPR151 polyclonal antibody

Catalog Number: PAB16377

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of GPR151.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to human GPR151.

Host: Rabbit

Reactivity: Gorilla, Human, Monkey

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: 3rd extracellular domain of human.

Form: Liquid

Purification: Immunoaffinity purification

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (20 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 134391

Gene Symbol: GPR151

Gene Alias: GALR4, GALRL, GPCR, PGR7