

Datasheet

FZD9 polyclonal antibody

Catalog Number: PAB16383

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of FZD9.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to human FZD9.

Host: Rabbit

Reactivity: Bovine,Dog,Human,Mouse,Rabbit,Rat

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: N-terminal extracellular domain of human.

Form: Liquid

Purification: Immunoaffinity purification

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (10 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8326

Gene Symbol: FZD9

Gene Alias: CD349, FZD3

Gene Summary: Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are

receptors for Wnt signaling proteins. The FZD9 gene is located within the Williams syndrome common deletion region of chromosome 7, and heterozygous deletion of the FZD9 gene may contribute to the Williams syndrome phenotype. FZD9 is expressed predominantly in brain, testis, eye, skeletal muscle, and kidney. [provided by RefSeq]