

Datasheet

HTR1E polyclonal antibody

Catalog Number: PAB16399

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of HTR1E.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to human HTR1E.

Host: Rabbit

Reactivity: Chimpanzee, Gorilla, Guinea pig, Horse, Human, Monkey

Applications: IHC-P
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: 3rd Cytoplasmic domain of human.

Form: Liquid

Purification: Immunoaffinity purification

Recommend Usage: Immunohistochemistry
(Formalin/PFA-fixed paraffin-embedded sections) (16-24 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 3354

Gene Symbol: HTR1E

Gene Alias: 5-HT1E