

Datasheet

PDE1C polyclonal antibody

Catalog Number: PAB16437

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of PDE1C.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to human PDE1C.

Host: Rabbit

Reactivity:

Bovine, Dog, Gorilla, Hamster, Horse, Human, Monkey, Mouse, Orangutan, Rabbit, Rat

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: C-terminal domain of human.

Form: Liquid

Purification: Immunoaffinity purification

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5137

Gene Symbol: PDE1C

Gene Alias: Hcam3