

Datasheet

ZNF397 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00084307-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human ZNF397 protein.

Immunogen: ZNF397 (NP_115723.1, 1 a.a. ~ 275 a.a) full-length human protein.

Sequence:

MAVESGVISTLIPQDPPEQELILVKVEDNFSWDEKFKQ
NGSTQSCQELFRQQFRKFCYQETPGPREALSRLQEL
CYQWLMPHELHTKEQILELLVLEQFLSILPEELQIWVQQ
HNPESGEEAVTLLEDLEREFDDPGQQVPASPQGPAV
PWKDLTCLRASQESTDIHLQPLKTQLKSWKPCLSPKS
DCENSETATKEGISEEKSQGLPQEPSFRGIKLSRPPKA
SSAIRWECVSPGSFPGDIAAEATHSTISCFAINLPTATI
LPSKNVNRKYFS

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84307

Gene Symbol: ZNF397

Gene Alias: MGC13250, ZNF47, ZSCAN15