

Datasheet

NTNG2 monoclonal antibody (M01), clone 4F11

Catalog Number: H00084628-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant NTNG2.

Clone Name: 4F11

Immunogen: NTNG2 (NP_115925, 21 a.a. ~ 119 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

ICKSWVTTDEGPTWEFYACQPKVMRLKDYVKVKVEPS
GITCGDPPERFCSHENPYLCSNECDASNPDLAHPRL
MFDKEEEGLATYWQSITWSRYPSP

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 84628

Gene Symbol: NTNG2

Gene Alias: KIAA0625, KIAA1857, LHLL9381, Lmnt2,
MGC21884, NTNG1, bA479K20.1