

Datasheet

ABHD1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00084696-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ABHD1 protein.

Immunogen: ABHD1 (AAH39576.2, 1 a.a. ~ 405 a.a) full-length human protein.

Sequence:

MLSSFLSPQNGTWADTFSLLLALAVALYLGYYWACVL
QRPRLVAGPQFLAFLEPHCSITTETFYPTLWCFEGRQLQ
SIFQVLLQSQPLVLYQSDILQTPDGGQLLLDWAKQPDS
SQDPDPTTQPIVLLPGITGSSQETVYLHLVNQALRDG
YQAVVFNNRGCGRGEELRTHRAFCASNTEDLETVVNHI
KHRYPQAPLLAVGISFGGILVLNHLAQARQAAGLVAAL
TLSACWDSFETTRSLETPLNSLLFNQPLTAGLCQLVER
NRKVIEKVVDIDFVLQARTIRQFDERYTSVAFGYQDCV
TYYKAASPRTKIDAIRIPVLYLSAADDPFSPVCALPIQAA
QHSPYVALLITARGGHIGFLEGLLPWQHWYMSRLLHQ
YAKAIFQDPEGLPDLRALLPSEDRNS

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84696

Gene Symbol: ABHD1

Gene Alias: FLJ36128, LABH1

Gene Summary: This gene is a member of the AB

hydrolase superfamily and encodes a protein with an alpha/beta hydrolase fold. This domain is common to a number of hydrolytic enzymes of widely differing phylogenetic origins and catalytic functions. [provided by RefSeq]