

Datasheet

SPSB2 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00084727-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human SPSB2 protein.

Immunogen: SPSB2 (NP_116030.1, 1 a.a. ~ 263 a.a) full-length human protein.

Sequence:

MGQTALAGGSSSTPTPQALYDLSCEGLEELLSAPP
PDLGAQRRHGWNPKDCSENIEVKEGGLYFERRPVAQ
STDGARGKRGYSRGLHAWAISWPLEQQRGTHAVVGVA
TALAPLQTDHYAALLGSNSESWGWDIGRGKLYHQSK
GPGAPQYPAGTQGEQLEVPERLLVLDMEEGTLGYAI
GGTYLGPAFRGLKGRTLYPAVSAVWGQCQVRIRYL
ERRAEPHSLHLSRLCVRHNLGDTRLGQVSALPLPPA
MKRYLLYQ

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84727

Gene Symbol: SPSB2

Gene Alias: GRCC9, MGC2519, SSB-2, SSB2

Gene Summary: This gene encodes a suppressor of cytokine signaling (SOCS) family member, and it belongs to the subfamily of proteins containing a central SPRY (repeats in splA and RyR) domain and a

C-terminal SOCS box. This gene is present in a gene-rich cluster on chromosome 12p13 in the vicinity of the CD4 antigen and triosephosphate isomerase genes. Alternative splicing results in multiple transcript variants. [provided by RefSeq]