

Datasheet

TUBA1C purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00084790-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human TUBA1C protein.

Immunogen: TUBA1C (NP_116093.1, 1 a.a. ~ 449 a.a) full-length human protein.

Sequence:

MRECISIHVGQAGVQIGNACWELYCLEHGIQPDGQMP
SDKTIGGGDDSFNTFFSETGAGKHVPRAVFVDLEPTVI
DEVRTGTYRQLFHPEQLITGKEDAANNYARGHYTIGK
EIIDLVLDRIKLADQCTGLQGFLVFHSFGGGTSGSFT
SLLMERLSVDYGGKSKLEFSIYPAPQVSTAVVEPYNSIL
TTHTTLEHSDCAFMVDNEAIYDICRRNLDIERPTYTNLN
RLISQIVSSITASLRFDGALNVDLTFEQTNLVPYPRIHFP
LATYAPVISAEKAYHEQLTVAEITNACFEPANQMVKCD
PRHGKYMCCLLYRGDVVPKDVNAAIATIKTKRTIQFV
DWCPTGFKVGINYQPPTVVPGGDLAKVQRAVCMLSN
TTAVAEAWARLDHKFDLMYAKRAFWHWYVGEGMEEG
EFSEAREDMAALEKDYEEVGADSADGEDEGEY

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84790

Gene Symbol: TUBA1C

Gene Alias: MGC10851, MGC14580, TUBA6, bcm948