

Datasheet

ZNF496 purified MaxPab mouse polyclonal antibody (B01P)

Gene Symbol: ZNF496

Gene Alias: MGC15548, NIZP1, ZKSCAN17

Catalog Number: H00084838-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF496 protein.

Immunogen: ZNF496 (NP_116141.1, 1 a.a. ~ 587 a.a) full-length human protein.

Sequence:

MPTALCPRVLAPKESEEPRKMRSPPGENPSPQGELP
SPESSRRLLFRFRYQEAAGPREALQRLWDLGGWLR
PERHTKEQILELLVLEQLAILPREIQSWVRAQEPESGE
QAVAAVEALEREPGRPWQWLKHCEDPVVIDDGDSPL
DQEQEQLPVEPHSDLAQNQDAQPITLAQCLGLPSRPP
SQLSGDPVLQDAFLLQEENVVRTQQTTLQLPPSRVS
PFKDMILCFSEEDWSLLDPAQTGFYGEFIIGEDYGVSM
PPNDLAAQPDLSQGEENEPRVPELQDLQGKEVPQVS
YLDSPSLQPFQVEERRKREELQVPEFQACPQTVVPQN
TYPAGGNPRSLNSLDEEVTIEIVLSSSGDEDSQHGPY
CTEELGSPTEKQRSLPASHRSSTEAGGEVQTSKKSIV
CPNCGKIFRWRVNFIRHLRSRREQKPHECSVCGELF
SDSEDLDGHLESHEAQKPYRCGACGKSFRLNSHLLSH
RRIHLQPDRLQPVEKREQAASEDADKGPKLENGKA
KLSFQCCECGKAFQRHDHLARHRSHFHLKDKARPFQ
CRYCVKSFTQNYDLLRHERLHMKRRSKQALNSY

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84838