

Datasheet

ZNF496 monoclonal antibody (M02), clone 4B1

Catalog Number: H00084838-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant ZNF496.

Clone Name: 4B1

Immunogen: ZNF496 (NP_116141, 485 a.a. ~ 585 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

HLQPDRLQPVEKREQAASEDADKGPKEPLENGKAKLS
FQCCECGKAFQRHDHLARHRSHFHLKDKARPFQCRY
CVKSFTQNYDLLRHERLHMKRRSKQALN

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 84838

Gene Symbol: ZNF496

Gene Alias: MGC15548, NIZP1, ZKSCAN17