

Datasheet

MINA MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00084864-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MINA protein.

Immunogen: MINA (AAH14928.1, 1 a.a. ~ 465 a.a)
full-length human protein.

Sequence:

MPKKAKPTGSGKEEGPAPCKQMKLEAAGGPSALNFD
SPSSLFESLISPIKTETFFKEFWQKPLLIQRDDPALAT
YYGSLFKLTDLKSLCSRGMYGRDVNVCRCVNGKKK
VLNKGDKAHFLQLRKDFDQKRATIQFHQPQRFKDELW
RIQEKLECYFGSLVGSNVYITPAGSQGLPPHYDDVEVF
ILQLEGEKHWRLYHPTVPLAREYSVEAEERIGRPVHEF
MLKPGDLLYFPRGTIHQADTPAGLAHSTHVTISTYQNN
SWGDFLLDTISGLVFD TAKEDVELRTGIPRQLLQVES
TTVATRRLSGFLRTLADRLEGTKELLSSDMKKDFIMHR
LPPYSAGDGAELSTPGGKLPRLDSVVRQLQFDHIVLTV
LPDQDQSDQETQEKMVYIYHSLKNSRETHMMGNEEET
EFHGLRFPLSHLDALKQIWNSPAISVKDLKLTDEEKE
SLVLSLWTECLIQVV

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84864

Gene Symbol: MINA

Gene Alias: DKFZp762O1912, FLJ14393, MDIG,

MINA53, NO52

Gene Summary: MINA is a c-Myc (MYC; MIM 190080) target gene that may play a role in cell proliferation or regulation of cell growth. (Tsuneoka et al., 2002 [PubMed 12091391]; Zhang et al., 2005 [PubMed 15897898]).[supplied by OMIM]