

Datasheet

C3orf39 purified MaxPab mouse polyclonal antibody (B01P)

Gene Symbol: C3orf39

Gene Alias: AGO61, FLJ14566

Catalog Number: H00084892-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human C3orf39 protein.

Immunogen: C3orf39 (NP_116195.2, 1 a.a. ~ 580 a.a) full-length human protein.

Sequence:

MHLSAVFNALLVSVLA AVLWKHVRLREHAATLEEEELAL
SRQATEPAPALRIDYPKALQILMEGGTHMVCTGRHTHT
DRICRFKWL CYSNEAE E FIFFHGNTSVMLPNLGSRRF
QPALLDLSTVEDHNTQYFNFVELPAAALRFMPKPVFVP
DVALIANRFNPDNLMHVFHDDLPLFYTLRQFPGLAHE
ARLFFMEGWGEGA HFDLYKLLSPKQPLLRAQLKTLGR
LLCFSHAFVGLSKITTWYQYGFVQPQGPKANILVSGNE
IRQFARFMTEKLNVSHTGVPLGEEYILVFSRTQNRLILN
EAELLALAQEFQMKT VTVSLEDHTFADVRLVSNAS
MLVSMHGAQLVTT LFLPRGATVVELFPYAVNPDHYTP
YKTLAMLPGMDLQYVAWRNMMPENTVTHPERPWDQ
GGITHLDRAEQARILQSREVPRHLCCRNPEWLFRIYQD
TKVDIPSLIQTIRRVVKGRPGPRKQKWTVGLYPGKVRE
ARCQASVHGASEARLTVSWQIPWNLKYLKVREVKYEV
WLQEQGENTYVPYILALQNHTFTENIKPFTTYLVWVRC
IFNKILLGPFADVLVCNT

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84892