

Datasheet

TNS4 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00084951-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human TNS4 protein.

Immunogen: TNS4 (AAH18706.1, 1 a.a. ~ 480 a.a) full-length human protein.

Sequence:

MGSKASSPHGLGSPLVASPRLEKRLGGLAPQRGSRIS
VLSASPVSDVSYMFGSSQSLHSSNSSHQSSRSLES
PANSSSLHSLGSVSLCTRPSTDFQAPRNPTLTMGQPR
TPHSPPLAKEHASSCPPSITNSMVDIPIVLINGCPEPGS
SPPQRTPGHQNSVQPGAASPSNPCPATRSNSQTLSD
APFTTCPEGPARDMQPTMKFVMDTSKYWFKPNITRE
QAIELLRKEEPGAFVIRDSSSYRGSFGLALKVQEVPA
AQSRPGEDSNDLIRHFLIESSAKGVHLKGADDEEYFGS
LSAFVCQHSIMALALPCKLTIPQRELGGADGASDSTDS
PASCQKKSAGCHTLYLSSVSVETLTGALAVQKAISTTF
ERDILPTPTVVHFKVTEQGITLTDVQRKVFFRRHYPLTT
LRFCGMDPEQRKWQKYCKPSWIFGFVAKSQTEPQEN
VCHLFAEYDMVQPASQVIGLVTALLQDAERM

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84951

Gene Symbol: TNS4

Gene Alias: CTEN, FLJ14950, PP14434