

Datasheet

MGC15875 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00085007-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MGC15875 protein.

Immunogen: MGC15875 (NP_116310.1, 1 a.a. ~ 175 a.a) full-length human protein.

Sequence:

MGKSIGNGHPVACVAATQPVARAFEATGVEYFNTFGG
SPVSCAVGLAVLNVLKEQLQDHATSVGSFLMQLLGQ
QKIKHPIVGDVRGVGLFIGVDLIKDEATRTPATEEAAYL
VSRLKENYVLLSTDGPGRNLIKFKPPMCFSLDNARQV
VAKLDAILTDMEEKVRSCETLRLQP

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 85007

Gene Symbol: AGXT2L2

Gene Alias: MGC117348, MGC15875, MGC45484