

Datasheet

PPIL4 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00085313-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human PPIL4 protein.

Immunogen: PPIL4 (NP_624311.1, 1 a.a. ~ 492 a.a) full-length human protein.

Sequence:

MAVLLETTGLGDVVIDLYTEERPRACLNFLKLCIKIYYNY
CLIHNVQRDFIIQTGDPTGTGRGGESIFGQLYGDQASF
FEAEKVPRIKHKKKGTVMVNNGSDQHGSQFLITTGE
NLDYLDGVHTVFGVTEGMDIHKINETFVDKDFVPYQ
DIRINHVTILDDPFDDPPDLLIPDRSPEPTREQLDSGRI
GADEEIDDFKGRSAEEVEEIKAEKEAKTQAILLEMVGD
LPDADIKPPENVLFVCKLNPVTTDEDLEIIFSRFGPIRSC
EVIRDWKTGESLCYAFIEFEKEEDCEKAFFKMDNVLID
DRRIHVDFSQSVAKVKWKGKGGKYTKSDFKEYEKEQ
DKPPNLVLKDKVKPKQDTKYDLILDEQAEDSKSSHSHT
SKKHKKKTHHCSEEKEDEDYMPIKNTNQDIYREMGFG
HYEEEEESCWEKQKSEKRDRTQNRSRSRERDGHY
SNSHKSKYQTDLYERERSKKRDRSRSPKSKDKEKSK
YR

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 85313

Gene Symbol: PPIL4

Gene Alias: HDCME13P

Gene Summary: This gene is a member of the cyclophilin family of peptidylprolyl isomerases. The cyclophilins are a highly conserved family, members of which play an important role in protein folding, immunosuppression by cyclosporin A, and infection of HIV-1 virions. [provided by RefSeq]