

Datasheet

RAGE monoclonal antibody, clone 1C5

Catalog Number: MAB7936

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against recombinant RAGE.

Clone Name: 1C5

Immunogen: Recombinant protein corresponding to human RAGE.

Host: Mouse

Reactivity: Human

Applications: IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Concentration: 250 ug/mL

Isotype: IgG, kappa

Recommend Usage: ELISA (0.01 ug/mL)

Flow Cytometry (1 ug/mL)

ICC (1 ug/mL)

Western Blot (0.1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% proclin, 2% Block Ace)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5891

Gene Symbol: RAGE

Gene Alias: MOK, RAGE1