

Datasheet

***Helicobacter pylori* CagA monoclonal antibody, clone HP-1811**

Catalog Number: MAB8529

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against partial recombinant *Helicobacter pylori*
CagA.

Clone Name: HP-1811

Immunogen: Recombinant protein corresponding to
C-terminus of *Helicobacter pylori* CagA.

Host: Mouse

Reactivity: *Helicobacter pylori*

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: MAbs HP-1811 recognizes epitope localized
between Asp562 and Gln795 in the GenBank FJ428215
amino acid sequence.

Form: Liquid

Isotype: IgG

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.1% sodium azide)

Storage Instruction: Store in the dark at 4°C. Do not
freeze.