

## Datasheet

### Influenza B virus monoclonal antibody, clone 15-12

**Catalog Number:** MAB8538

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against influenza B virus.

**Clone Name:** 15-12

**Immunogen:** Native purified influenza B virus.

**Host:** Mouse

**Reactivity:** Viruses

**Applications:** ELISA

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Nucleoprotein of influenza virus type B.

**Form:** Liquid

**Isotype:** IgG2a

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (0.1% sodium azide, 0.1% gelatin)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.