

Datasheet

C17orf57 MaxPab mouse polyclonal antibody (B01)

Gene Alias: FLJ40342

Catalog Number: H00124989-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human C17orf57 protein.

Immunogen: C17orf57 (AAH37876, 1 a.a. ~ 531 a.a) full-length human protein.

Sequence:

MENDDLESKRPKNTWQIRKFLGGVGSSNVGVQEPYS
KNGINFKKHSEKGEIHDSKSKPQSLKSSTLSKSLDKS
DISSIPKLQKPAVRKHSSLQKQVSSTEKTAISTLENFCE
AISKLQENYIAAEELQSILPSTGINLLDEEFQKIVTDTSR
NENGMVELDDFVNALAKERSFPECNALPGVIKAIKIK
DKNVDYEDLNTCLQNFGIYLSKPEFKKITELTEAGETKK
VNFKEFIDTMMSNTECFSEKLVLPDAIETLDDLKMETM
SVSDLWNTLSSLNSNLKKDEFLLAEELVTVDGDKVQ
FEFAKVVVRNMRDAARLEELQEVVLAADLLEGDMIAG
KNLEDFLRNVGIKSPKEEVEKILQSDFVSEDNMVNIKD
CMRALRDTQKFSNYIDFRKEASNLKLPKVNEIKEAANIL
SHVDNGKIGIPDLEHALKCLNVNLTEEDFNEALNCCNV
SDNMEVDLKDFLMKMKESPHFQKSKATQILLATTQILQ
NDLVDVSDLKTLMDKDLHTANAILTVMLRHVPEHG

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 124989

Gene Symbol: C17orf57