

Datasheet

MGC33839 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00125875-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MGC33839 protein.

Immunogen: MGC33839 (NP_689566, 1 a.a. ~ 167 a.a) full-length human protein.

Sequence:

MGVKRSLQSGGILLSLVANVLMVLSTATNYWTRQQEG
HSLGWQECNHGICSSIPCQTTLAVTVACMVLAVGVGV
VGMVMGLRIRCDEGESLRGQTTS AFLFLGLLLLLTALI
GYTVKNAWKNNVFFSWSYFSGWLALPFSILAGFCFL
ADMIMQSTDAISGFPVCL

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 125875

Gene Symbol: CLDND2

Gene Alias: MGC33839