

Datasheet

COX6B2 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00125965-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human COX6B2 protein.

Immunogen: COX6B2 (NP_653214.2, 1 a.a. ~ 88 a.a) full-length human protein.

Sequence:

MLDVEAQEPPKGGKSTPPFDPRFPSQNGIRNCYQNF
LDYHRCLKTRTRRGKSTQPCEYYFRVYHSLCPISWVE
SWNEQIKNGIFAGKI

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 125965

Gene Symbol: COX6B2

Gene Alias: COXVIB2, MGC119094