

Datasheet

MGC29891 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00126626-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MGC29891 protein.

Immunogen: MGC29891 (NP_653219.1, 1 a.a. ~ 448 a.a) full-length human protein.

Sequence:

MSLVDLGKRLLEAARKGQDDEVRTL MANGAPFTTDW
LGTSPHLAAQYGHYSTAEVLLRAGVSRDARTKVDRT
PLHMAAADGHAHIVELLVRNGADVNAKMLKMTALH
WATERHHRDVVELLIKYGADVHAFSKFDKSAFDIALEK
NNAEILVILQEAMQNQVNVNPERANPVTDPVSMAAPFI
FTSGEVVNLASLISSTNTKTTSGDPHASTVQFSNSTTS
VLATLAALAEASVPLSN SHRATANTEEIIIEGNSVDSSIQ
QVMGSGGQRVITIVTDGVPLGNIQTSIPTGGIGQPFIVT
VQDGGQVLTVPAGKVAEETVIKEEEEEKLPLTKKPRIG
EKTNSVEESKEGNERELLQQQLQEANRRAEYRHQL
LKKEQAEQYRLKLEAIARQQPNGVDFTMVEEVAEVD
AVVVTGELEERETKVTGSAGTTEPHTRVSMATVSS

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 126626

Gene Symbol: GABPB2

Gene Alias: MGC29891