

Datasheet

LRRC44 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00127255-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LRRC44 protein.

Immunogen: LRRC44 (AAH20778, 1 a.a. ~ 199 a.a) full-length human protein.

Sequence:

MFHGTVTEELTSHEEWSHYNENIREGQKDFVFKFNG
LHLKSMENLQSCISLRVCIFSNNFITDIHPLQSCIKLIKLD
LHGNQIKSLPNTKFWNGLKNLKLlyLHDNGFAKLKNIC
VLSACPTLIALTMFDCPVSLKKGyRHVLVNSIWPLKAL
DHHVISDEEIIQNWHLPERFKACNHRLFFNFCPALRKE
KMKHSEV

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 127255

Gene Symbol: LRR1Q3

Gene Alias: LRRC44, MGC22773