

Datasheet

UHMK1 polyclonal antibody (A01)

Catalog Number: H00127933-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant UHMK1.

Immunogen: UHMK1 (AAH26046, 1 a.a. ~ 105 a.a)
partial recombinant protein with GST tag.

Sequence:

MAGSGCAWGAEPFRFLEAFGRWLQVQSRLGSGSSA
SVYRVRCCGNPGSPPGALKQFLPPGTTGAAASAAEY
GFRKERAALQLQGHNRNIVTLYGVFTIHFSPNVP

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 127933

Gene Symbol: UHMK1

Gene Alias: KIS, Kist