

Datasheet

SIRPD purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00128646-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human SIRPD protein.

Immunogen: SIRPD (NP_848555.1, 1 a.a. ~ 197 a.a) full-length human protein.

Sequence:

MPIPASPLHPPLPSLLLYLLELAGVTHVFHVQQTEMS
QTVSTGESIILSCSVPTLPNGPVLWFKGTGPNRKLIY
NFKQGNFPRVKEIGDTTKPGNTDFSTRIREISLADAGT
YYCVKFIKGRAIKEYQSGRGTQVFVTEQNPRPPKNRP
AGRAGSRAHDAHTCLSALPERNSTNYFVQPCCCLRL
LGLTGLLSK

Host: Mouse

Reactivity: Human

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 128646

Gene Symbol: SIRPD

Gene Alias: PTPNS1L2, dJ576H24.4