

Datasheet

LOC129530 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00129530-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LOC129530 protein.

Immunogen: LOC129530 (NP_777558.1, 1 a.a. ~ 194 a.a) full-length human protein.

Sequence:

MSALWLLLGLLALMDLSESSNWGCYGNISLDTPGAS
CGIGRRHGLNYCGVRASERLAEIDMPYLLKYQPMMQT
IGQKYCMDPAVIAGVLSRKSPGDKILVNMGDRTSMVQ
DPGSQAPTWISESQVSQTTEVLTTRIKEIQRRFPTWT
PDQYLRGGLCAYSGGAGYVRSSQDLSCDFCNDVLAR
AKYLKRHGF

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 129530

Gene Symbol: LYG1

Gene Alias: SALW1939