

Datasheet

KCTD18 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00130535-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human KCTD18 protein.

Immunogen: KCTD18 (NP_689600.2, 1 a.a. ~ 426 a.a) full-length human protein.

Sequence:

MEGHKAEVEVLDVLRNLNVGGCIYTARRESLCRFKDSM
LASMFSGRFPLKTDESGACVIDRDGRLFKYLLDYLGHE
VQIPTDEQTRIALQEEADYFGIPYPYSLSDHLANEMET
YSLRSNIELKKALTDFCDSYGLVCNKPTVWVLHYLNTS
GASCESRIIGVYATKTDGTDIEKQLGGRIHSGIFKRE
AGNNVQYIWSYYSVAELKKMMDAFDAWEGKGVSYW
RVPHELIECWTLERPLLGLSRHMAPIRKRRLLITFNEAD
ESVNYKTGPKPVRFLGPSTSTQIKVKNSASVTVSPASA
IQTSAGATANRFQSGSRRKAAQRSAPSRATALVGTGA
PGHPQASPGAASAENGTHLPPAKVLLSDKKPTPQRV
IKLKRTPLCATAPCLPSTATRQANSLKPLPGEAARAL
GVRTENGKNKGN

Host: Mouse

Reactivity: Human

Applications: IF, WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 130535

Gene Symbol: KCTD18

Gene Alias: 6530404F10Rik, FLJ31322, FLJ37818