

Datasheet

GPR175 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00131601-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human GPR175 protein.

Immunogen: GPR175 (AAH50711.1, 1 a.a. ~ 373 a.a) full-length human protein.

Sequence:

MDTLEEVWANGSTALPPPLAPNISVPHRCLLLLYEDI
GTSRVRYWDLILLIPNVLFLLWKLPSARAKIRITSSPI
FITFYILVFVVALVGIAVVSMVTSTNAATVADKILWE
ITRFFLLAIELSVIILGLAFGHLESKSSIKRVLAITTVLSLA
YSVTQGTLEILYPDAHLSAEDFNIYGHGGRQFWLVSS
CFFFLVYSLVILPKTPLKERISLPSRRSFYVYAGILALL
NLLQGLGSVLLCFDIEGLCCVDATTFLYFSFFAPLIYVA
FLRGFFGSEPKILFSYKCQVDETEEPDVHLPQPYAVAR
REGLEAAGAAGASAASYSSTQFDSAGGVAYLDDIASM
PCHTGSINSTDSERWKAINA

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 131601

Gene Symbol: GPR175

Gene Alias: FLJ32197, TPRA40