

Datasheet

FLJ39155 polyclonal antibody (A01)

Catalog Number: H00133584-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant FLJ39155.

Immunogen: FLJ39155 (NP_689616, 914 a.a. ~ 1001 a.a) partial recombinant protein with GST tag.

Sequence:

DGRWHRVKAVRDGQSGKITVDDYGARTGKSPGMMR
QLNINGALYVGGMKEIALHTNRQYMRGLVGCISHFTLS
TDYHISLVEDAVDGK

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 133584

Gene Symbol: EGFLAM

Gene Alias: AGRINL, AGRNL, FLJ39155