

Datasheet

JMY monoclonal antibody (M01), clone 2D7

Catalog Number: H00133746-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant JMY.

Clone Name: 2D7

Immunogen: JMY (NP_689618, 535 a.a. ~ 634 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

SPMDEVLASLKRGSFHLKKVEQRTLPPFPDEDDSNIL
AQIRKGVKLKKVQKVDLRESFTLLPDTDPLTRSIHEALR
RIKEASPESEDEEEALPCTDWEN

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 133746

Gene Symbol: JMY

Gene Alias: FLJ37870, MGC163496