

Datasheet

GST tag monoclonal antibody, clone GST.B6/G2R

Catalog Number: MAB9750

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against GST.

Clone Name: GST.B6/G2R

Immunogen: GST protein.

Host: Mouse

Applications: Dot, ELISA, IHC, IP, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: Recognizes native and denatured forms of
purified GST or GST fusion proteins.

Form: Liquid

Purification: Protein A purification

Concentration: 1 mg/mL

Isotype: IgG2a

Recommend Usage: Western Blot (1:1000-1:3000)
Immunohistochemistry (1:500-1:2000)
The optimal working dilution should be determined by
the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (0.05% sodium
azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.