

Datasheet

Itgam monoclonal antibody, clone M1/70 (APC/Cy7)

Catalog Number: MAB9868

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against Itgam.

Clone Name: M1/70

Immunogen: Native Itgam from mouse splenic T cells and splenocytes.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: APC/Cy7

Isotype: IgG2b, kappa

Recommend Usage: Flow Cytometry (? 0.3 ug/10⁶ cells)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 16409

Gene Symbol: Itgam

Gene Alias: CD11b/CD18, CR3, CR3A, Cd11b, F730045J24Rik, Ly-40, MAC1, Mac-1, Mac-1a