

## Datasheet

### OCIAD1 polyclonal antibody

**Catalog Number:** PAB16855

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of OCIAD1.

**Immunogen:** A synthetic peptide corresponding to C-terminus 16 amino acids of human OCIAD1.

**Host:** Rabbit

**Reactivity:** Human

**Applications:** ELISA, WB-Ce  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** Western Blot (1-2 ug/mL)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.02% sodium azide)

**Storage Instruction:** Store at 4°C for three months. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 54940

**Gene Symbol:** OCIAD1

**Gene Alias:** Asrij, FLJ20455, MGC111072, OCIA, TPA018