

Datasheet

OCIAD2 polyclonal antibody

Catalog Number: PAB16857

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of OCIAD2.

Immunogen: A synthetic peptide corresponding to N-terminus 17 amino acids of human OCIAD2.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: ELISA, IF, WB-Ce
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (0.5-1 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 132299

Gene Symbol: OCIAD2

Gene Alias: DKFZp686C03164, MGC45416