

Datasheet

SRPX2 polyclonal antibody

Catalog Number: PAB16859

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of SRPX2.

Immunogen: A synthetic peptide corresponding to N-terminus 16 amino acids of human SRPX2.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: ELISA, IHC, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (1-2 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 27286

Gene Symbol: SRPX2

Gene Alias: BPP, CBPS, RESDX, SRPUL