

## Datasheet

### Trpv3 polyclonal antibody

**Catalog Number:** PAB16890

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Trpv3.

**Immunogen:** A synthetic peptide (conjugated with KLH) corresponding to Trpv3.

**Sequence:** YRPREEEAI

**Host:** Rabbit

**Reactivity:** Mouse

**Applications:** EIA, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Lyophilized

**Recommend Usage:** EIA (1:5000)

Western Blot (1:2000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 10 mM PBS

**Storage Instruction:** Store at -20°C on dry atmosphere.

After reconstitution with deionized water, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 246788

**Gene Symbol:** Trpv3

**Gene Alias:** 1110036l10Rik, AI644701, MGC124324, MGC124325, Nh, VRL3