

Datasheet

BRP44 polyclonal antibody

Catalog Number: PAB17634

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of BRP44.

Immunogen: A synthetic peptide corresponding to internal of human BRP44.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: ELISA, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody detects endogenous levels of total BRP44 protein.

Form: Liquid

Recommend Usage: Immunohistochemistry

(1:50-1:100)

ELISA (1:40000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 25874

Gene Symbol: BRP44

Gene Alias: DKFZp564B167, MGC125752, MGC125753